


Quick Reference Guide



PURITY - RECOVERY - THROUGHPUT TRILEMMA



FRACTION COLLECTION: A TRADE-OFF BETWEEN RECOVERY AND PURITY


[illegible]

Figure 1 consists of two panels. The left panel is a schematic of a fly on a stick. The fly's body is divided into three regions: the head (ground out), the thorax (base hit), and the abdomen (pop fly). Above the fly, the corresponding K_0 ranges are indicated: $K_0 < 0.4$ for ground out, $0.4 < K_0 < 2.5$ for base hit, and $K_0 > 2.5$ for pop fly. The right panel is a graph of the detector signal versus elution volume. The signal shows several peaks. A green oval highlights the region where the signal is highest, corresponding to the base hit region. Arrows point from the K_0 values 0.4, 1, and 2.5 to the corresponding peaks in the signal.

TERNARY DIAGRAMS




Diagram illustrating the relative hydrophobicity of various water-soluble polymers, categorized by their position on a scale from Lipophilic (left) to Hydrophilic (right).


The scale is represented by a horizontal bar with a green arrow pointing left (Lipophilic) and a green arrow pointing right (Hydrophilic).

The polymers are arranged in a staggered fashion, indicating their relative positions on the scale:

- HEMWat (top, dark blue bar)
- HterAcWat (light blue bar)
- EBUWAT (light blue bar)
- terAcWAT (light blue bar)
- ChMWAT (light blue bar)
- ARIZONA (HdEMWAT) (bottom, dark blue bar)

Arizona
HpEMW

Construct Fractogram



Selectivity	Theoretical plate number (Efficiency)	Resolution
$\alpha = \frac{k_2}{k_1} = \frac{t_{R,2} - t_0}{t_{R,1} - t_0} \geq 1.05$	$N = 4 \left(\frac{V_R}{W_{0.68}} \right)^2 = 16 \left(\frac{V_R}{W_S} \right)^2$	$R_S = \frac{V_{R,2} - V_{R,1}}{W_{R,2} + W_{R,1}} \geq 1.5$
	$W_{0.68}$ peak width at 60% of peak height W_S peak width at baseline	$R_S = \frac{\sqrt{N}}{4} \frac{\alpha - 1}{\alpha + 1} \frac{k}{1 + k}$
		$R_S = S_F \frac{\sqrt{N}}{4} \frac{(K_{D,2} - K_{D,1})}{1 - S_F(1 - \frac{K_{D,1}}{2} - \frac{K_{D,2}}{2})}$
Pressure drop on the rotor		Centrifugal force
$P_d = n (d\omega^2 Rh + \nu r F)$		$g = R\omega^2$
Hydrostatic term Hydrodynamic term		
n is the number of cells		
d is the phase density difference		
ω is the rotor spin rate		
R is the rotor radius		
h is the stationary phase height in a cell		
ν is the mobile phase viscosity		
F is a constant for the geometry of the apparatus		
S_F is the mobile phase flow rate		

Hydrostatic term	Hydrodynamic term
n	is the number of cells
$\Delta\omega$	is the phase density difference
ω	is the rotor spin rate
R	is the rotor radius
h	is the stationary phase height in a cell
v	is the mobile phase viscosity
γ	is a constant for the geometry of the apparatus
F	is the mobile phase flow rate