

CPC & HPLC

A TECHNOLOGY LEVEL COMPARISON

One Industrial Challenge, Two Separation Logics



CHAPTER 1

Chromatography is often unavoidable in the industrial purification of complex, high-value molecules, especially when structurally similar impurities resist crystallization and simple extraction [1–3]. Industrial downstream purification is driven by the need for high selectivity, yet selectivity (α) is not an abstract property of the target alone; it is created by the physicochemical environment of the separation system [4–9]. At preparative scale, the decisive question is not simply whether a target can be resolved from its impurities, but whether the retention mechanism, transport regime, loading behavior and scale-up burden of the separation system remain compatible with throughput, recovery, productivity and process control. For that reason, a comparison should begin at the level of separation physics, rather than at the level of named equipment platforms [4–9].

Chromatography emerged in the early twentieth century and evolved from a laboratory method into a core industrial purification technology as demands for selectivity, reproducibility and scale increased. Modern high-pressure liquid chromatography (HPLC) marked a decisive step in that evolution. Csaba Horváth, a Hungarian-born scientist, contributed directly to that transition through the development of the first modern HPLC and through the demonstration that rapid, high-efficiency liquid separations could be achieved by combining high pressure with small-particle column packings [10, 11].

Complex active pharmaceutical ingredient (API) streams rarely contain a single dominant impurity class. Instead, the target is embedded in a chemically crowded matrix comprising closely related structural derivatives, hydrolysis products, deprotection byproducts, oligomeric forms [1–3]. Industrial purification therefore depends on selective, controllable, well-controlled chromatographic operations, frequently arranged as a multistep cascade, combining a first enrichment or burden-reduction step with a later polishing step [3].

1. Governing descriptors of retention behavior and transport constraints

1.1 Liquid-solid chromatographic environment

Any preparative chromatography system is ultimately judged by the same governing descriptors: retention factor (k'), resolution (R_s), efficiency (N , theoretical plate number) and selectivity factor (α). Retention factor (Eq. 1.) is given by

$$k' = \frac{t_R - t_0}{t_0} \quad (\text{Eq. 1.})$$

where t_R is the retention time of the solute and t_0 is the dead time (hold-up time) of an unretained species. Thus, k' expresses how strongly a solute is retained relative to the passage of the mobile phase through the system. A value close to zero indicates little retention, whereas larger values indicate stronger retention, longer residence time and therefore longer retention times in the separation environment; under otherwise fixed operating conditions, this often also implies a longer cycle time.

For more than one solute, chromatographic discrimination is commonly expressed through the selectivity factor (Eq. 2.)

$$\alpha = \frac{k'_2}{k'_1} \quad (k'_2 > k'_1) \quad (\text{Eq. 2.})$$

where k'_1 and k'_2 are the retention factors of the less retained and more retained component, respectively. When $\alpha = 1$, the two solutes are not chromatographically distinguished, as α increases above unity, separation becomes progressively feasible. This is why even modest changes in selectivity can have major practical consequences.

Chromatographic efficiency (Eq. 3. and 4.) is commonly expressed through the theoretical plate number (N), which describes how strongly solute bands broaden as they migrate through the system.

$$N = 16 \left(\frac{t_R}{W_b} \right)^2 \quad (\text{Eq. 3.})$$

when the baseline peak width, W_b , is used, or as

$$N = 5.54 \left(\frac{t_R}{W_{1/2}} \right)^2 \quad (\text{Eq. 4.})$$

when the peak width at half maximum, $W_{1/2}$, is used (Eq. 4.). In other words, efficiency describes how effectively a column or separation system limits dispersion while the solute band migrates.

The theoretical plate number is directly related to plate height (Eq. 5.), H , through

$$N = \frac{L}{H} \quad (\text{Eq. 5.})$$

where L is the effective separation length and H is the height equivalent to a theoretical plate. Plate height is therefore a measure of how much band broadening occurs per unit length of the separation path. Smaller H values correspond to more efficient separations, whereas larger H values indicate greater dispersion and poorer efficiency. A higher N corresponds to narrower peaks at a given retention and therefore to a greater ability to separate closely eluting components. In packed-bed systems, N is strongly influenced by particle size, longitudinal diffusion, eddy dispersion, mass-transfer resistance and is therefore closely linked to both column architecture and operating conditions.

In packed-bed chromatography, plate height is commonly interpreted through the *van Deemter* relationship (Eq. 6.),

$$H = A + \frac{B}{u} + Cu \quad (\text{Eq. 6.})$$

where u is the linear mobile phase velocity. The A -term describes eddy dispersion arising from multiple flow pathways through the packed bed, the B/u -term represents longitudinal diffusion and the Cu -term accounts for finite mass transfer kinetics between phases. This relationship makes clear that efficiency is not a purely geometric property, but a transport-dependent outcome shaped by bed structure, flow regime and phase-transfer resistance. For that reason, attempts to increase productivity by increasing flow rate do not improve efficiency indefinitely, beyond a certain point, mass transfer limitations dominate and plate height increases.

This has direct consequences for the engineering logic of preparative HPLC (high performance liquid chromatography). High efficiency is commonly pursued through small particles and tightly packed beds, but both choices increase hydraulic resistance and make pressure drop a central scale-up burden. As column diameter and throughput increase, the bed must remain mechanically uniform under substantially greater loading and flow demand, while still limiting void formation, channeling and asymmetry drift. This is precisely why dynamic axial compression (DAC) column technology became important in preparative HPLC, not because it changes the chromatographic mechanism, but because it helps preserve bed integrity and efficiency as packed-bed operation is pushed toward industrial scale.

The combined effect of retention, selectivity and efficiency on peak separation is commonly expressed by the approximate resolution equation (Eq. 7.)

$$R_s \approx \frac{\sqrt{N}}{4} \cdot \frac{\alpha - 1}{\alpha} \cdot \frac{k'}{1 + k'} \quad (\text{Eq. 7.})$$

where N represents theoretical plate number and k' is the representative retention factor of the solute pair. This expression makes clear that resolution is not controlled by a single property, but by three coupled contributions: column efficiency (N), chemical discrimination between the solutes (α) and the degree of retention (k'). In practical terms, difficult separations in HPLC are often improved by increasing efficiency, that is, by generating narrower peaks through higher plate numbers.

In one broad class of preparative systems, the solute migrates in a flowing liquid while repeatedly interacting with a porous solid stationary phase. Retention is therefore generated in a heterogeneous interfacial field defined not only by the mobile phase and bonded stationary phase chemistry, but also by residual surface functionalities, pore structure and hardware surfaces. The observed chromatographic behavior is not a pure equilibrium partitioning phenomenon, but a composite outcome of solvation, adsorption-desorption equilibria, electrostatic contributions, pore accessibility, local surface heterogeneity and mass-transfer resistance. This has direct consequences for selectivity: discrimination between closely related solutes is not governed solely by their intrinsic molecular properties, but by how differently those properties are expressed at the solid-liquid interface. Residual silanol groups, ligand density, pore topology, adsorbed solvent layers, local ionization phenomena and hardware interactions may therefore all contribute to effective retention and selectivity. In such systems, selectivity can be highly powerful, but it remains structurally coupled to the surface chemistry and transport properties of the packed bed.

Optimizing the performance of solid stationary phase has its limitation due to cumbersome control of particle morphology, porosity and surface chemistry.

1.2 Liquid-liquid environment

By contrast, in liquid-liquid preparative systems the retention problem is reformulated entirely: the solute is not retained by repeated interaction with a solid surface, but by preferential distribution between two immiscible liquid phases. The governing descriptor (Eq. 8.) is therefore the partition coefficient, K_D ,

$$K_D = \frac{C_S}{C_M} \quad (\text{Eq. 8.})$$

where C_S and C_M are the equilibrium concentrations of the solute in the stationary and mobile liquid phases, respectively. In this case, retention is governed primarily by the free-energy difference associated with mass transfer between the two liquid phases. The dominant contributions are therefore differential solvation, polarity, H-bonding balance, ionization state and specific intermolecular associations within one phase or the other.

Because this distribution is an equilibrium property, it is directly linked to the *Gibbs* free energy of transfer between the two phases (Eq. 9.):

$$\Delta G_{\text{transfer}} = -RT \ln K_D \quad (\text{Eq. 9.})$$

where $\Delta G_{\text{transfer}}$ is the *Gibbs* free energy associated with transfer of the solute from the mobile liquid phase to the stationary liquid phase, R is the universal gas constant, T is the absolute temperature and $\ln K_D$ is the natural logarithm of the partition coefficient. This relationship makes it explicit that liquid-liquid retention is thermodynamic in origin. A solute is retained not because it repeatedly interacts with a solid surface, but because mass transfer between two immiscible phases is associated with a well-defined free energy difference. Selectivity therefore reflects differences in *Gibbs* free energy of transfer between solutes, i.e. differences in how favorably each component is solvated by the two competing liquid environments.

However, thermodynamic preference alone is not sufficient to generate chromatographic retention. For the partition equilibrium to be expressed chromatographically, the solute must repeatedly transfer across the liquid-liquid interface while the two phases are brought into dynamic contact inside the rotor. The practical retention behavior of the system therefore depends not only on K_D , but also on the extent of interfacial mass transfer and on how much stationary liquid phase is retained in the column. In this sense, K_D defines the thermodynamic driving force of retention, whereas stationary phase retention (S_f) determines how strongly that retention can be realized under operating conditions.

Chromatographic retention in CPC (centrifugal partition chromatography) is therefore governed not only by partition thermodynamics, but also by how much stationary liquid phase is retained in the column (more specifically in the rotor unit). For a liquid–liquid chromatographic instrument, as CPC, the stationary phase retention (Eq. 10.), S_f is commonly expressed as

$$S_f = \frac{V_S}{V_C} \quad (\text{Eq. 10.})$$

where V_S is the stationary phase volume retained in the rotor or column and V_C is the total column volume. Since the mobile-phase hold-up is correspondingly reduced as V_S increases, S_f determines the effective phase ratio available for partitioning. Under idealized conditions, the retention volume (V_R , Eq. 11.) may be approximated as

$$V_R = V_M + K_D V_S \quad (\text{Eq. 11.})$$

where V_M is the mobile phase hold-up volume. Chromatographic retention in liquid–liquid systems therefore depends jointly on partition thermodynamics and retained stationary phase volume: K_D determines the preference of the solute for one phase over the other, while S_f determines how much stationary liquid is available for that preference to be expressed chromatographically.

This has direct implications for process design. Maximizing S_f is not a thermodynamic goal, but an operational one: a higher retained stationary phase fraction increases the effective phase ratio, strengthens retention for a given K_D , improves loading capacity and reduces the risk that stationary phase loss (column bleed) will erode separation performance. For that reason, stationary phase retention is a central engineering variable in CPC, shaped by rotor design, rotational speed, flow rate, viscosity, density difference and interfacial tension between the two phases.

This distinction is also structural. In packed-bed systems, part of the column volume is occupied by the solid support itself and is therefore unavailable for liquid–liquid partitioning. By contrast, in CPC the column volume is distributed only between the mobile and stationary liquid phases, with no inert solid skeleton consuming separation volume. As a result, the effective phase ratio is governed directly by stationary phase retention rather than by the fixed geometric penalty of a porous support (Fig. 1.). This is one reason why maximizing S_f is so important in CPC: it determines how much of the rotor volume is actually available for thermodynamically driven partitioning.




Figure 1. Schematic comparison of phase-volume distribution in packed-bed liquid chromatography and support-free liquid–liquid chromatography. In packed-bed systems, part of the column volume is occupied by solid support, reducing the volume available for stationary and mobile liquid phase interactions. In CPC, the column volume is shared only between the two liquid phases, so the effective phase ratio depends directly on stationary phase retention.

Unlike in packed-bed systems, retention behavior is therefore linked more directly to measurable phase partition and retained stationary phase fraction than to a complex superposition of surface-mediated physicochemical events. This difference has important implications for selectivity engineering.

In liquid–liquid systems, selectivity is generated mainly by solvent design. Composition, phase ratio, ionic strength, pH and dissolved selectors or ligands can shift the partition coefficient and thereby alter the relative retention of solutes. The central design variable is therefore not only column chemistry, but the thermodynamic and hydrodynamic environment in which mass transfer between the two immiscible phases occurs.

1.3. Comparison of the interactions on a molecular level

This contrast becomes industrially decisive when the same impurity problem is evaluated not only in terms of nominal resolution, but in terms of how retention, selectivity, efficiency and transport are physically generated.

The contrast between the separation environments is thus fundamental. In liquid-solid chromatography, retention and selectivity are generated in a structured interfacial field with significant transport constraints imposed by the packed bed. In liquid-liquid chromatography, retention and selectivity are generated by phase partitioning within a matrix-free environment, where the main design leverage lies in solvation thermodynamics and stationary phase retention. **The first environment tends to favor very high efficiency when bed structure and transport can be tightly controlled. The second offers a different route to performance by changing how selectivity itself is generated.** Although, solvent system development for LLC is as time consuming as for HPLC, once it is found, however, CPC scale-up is typically straightforward and close to linear, because the same partition logic, stationary phase retention and phase behavior can be carried across rotor scales.

In industrial practice, these two separation environments are embodied by the two major preparative platforms discussed: HPLC, which operates through a liquid-solid separation architecture and centrifugal partition chromatography (CPC), which operates through a liquid-liquid separation architecture.

The comparison between them should therefore begin not with instrument labels, but with the fact that they solve chromatographic problems through different separation logics, different transport regimes and different routes to selectivity.

As defined in Eq. 2, the selectivity factor α describes how differently two solutes are retained by the separation system. When α is close to unity, practical separation becomes difficult, as α increases, chromatographic discrimination becomes progressively more feasible. In liquid-liquid systems and under comparable phase-ratio conditions, α may be approximated by the ratio of partition coefficients of the two solutes. This makes clear why solvent system design is so critical in CPC: composition, phase ratio, pH, ionic strength, additives, dissolved selectors and temperature can all shift partition behavior and thereby alter selectivity. Proper CPC method selection is therefore not based on finding acceptable K_D values alone, but on finding a system that simultaneously provides a workable K_D range and sufficient selectivity for the compound of interest (Col) and its impurities that defines the purification problem.

2. Selectivity engines and accessible operation modes

Once the retention environment is defined, the next question is the accessible mode space. In liquid-solid preparative chromatography, the dominant interaction field is largely prescribed by stationary phase chemistry. The developer selects among hydrophobic, polar adsorption, ionic, zwitterionic, chiral, affinity, mixed-mode or exclusion mechanisms and then tunes gradient profile, pH, ionic strength, temperature and modifier composition within that structural framework. Before comparing mode space, however, one fundamental difference in separation capacity must be stated plainly. Preparative HPLC, particularly with modern sub-2- μm or superficially porous particles, routinely generates 10,000–100,000 theoretical plates under optimized conditions. CPC, by contrast, typically achieves 200–3,000 theoretical plates under standard operating conditions [7, 9, 22].




Figure 2. Theoretical plate number (Eq. 3.) required to achieve a fixed resolution ($R_s = 1.0$, solid; $R_s = 1.5$, dashed) as a function of selectivity, from the resolution equation (Eq. 4). The hyperbolic divergence as $\alpha \rightarrow 1$ shows that efficiency requirements become prohibitive for low-selectivity separations. Conversely, the curve flattens at higher selectivity: once α is sufficiently large, the required plate number falls to a few hundred and becomes largely insensitive to further gains in efficiency, so a high-selectivity separation can be resolved with a modest plate count, and additional plates yield diminishing returns.

Projected onto the resolution requirement (Fig. 2.), this gap maps directly onto selectivity. At a target resolution of $R_s = 1.0$, the 10^4 – 10^5 plates available to modern preparative HPLC suffice to baseline-resolve solute pairs with selectivities as low as $\alpha \approx 1.02$ – 1.05 , whereas the 200–3,000 plates typical of CPC require α to reach roughly 1.1–1.5 to achieve the same resolution. The hyperbolic divergence as $\alpha \rightarrow 1$ shows that efficiency requirements become prohibitive for low-selectivity separations; conversely, at higher selectivity the requirement collapses to a few hundred plates and flattens, so favorable selectivity can compensate for a low plate number.

The unresolved-peak problem is therefore answered differently by the two platforms. In preparative HPLC, the usual route is to increase efficiency, that is, to generate narrower peaks through higher plate number and kinetic optimization of the packed bed. In CPC, the usual route is to increase selectivity, that is, to change the relative partitioning of the solutes through solvent system design. In other words, preparative HPLC often solves difficult separations by pushing N upward, whereas CPC solves them by pushing α upward (Eq. 7.). This is why the lower theoretical plate number of CPC is not, by itself, a disqualifying limitation: if solvent system design produces sufficiently favorable selectivity, the separation problem can be solved by a different route (Fig. 3.).




Figure 3. Two routes to improved separation performance for the same unresolved impurity pair: preparative HPLC primarily by increasing efficiency (N) through kinetic optimization and CPC primarily by increasing selectivity (α) through solvent-system design.

This means that CPC cannot replicate HPLC performance in separations where selectivity is low (α close to 1.0) and resolution therefore depends critically on high N . When α is sufficiently large ($\alpha \geq 1.5$), or when the selectivity achievable by solvent system design exceeds what is accessible in the liquid–solid environment, the efficiency gap becomes practically irrelevant. This distinction between efficiency-dependent and selectivity-driven separations is the correct frame for platform selection and it will be developed further in the performance comparison that follows.

In liquid-liquid chromatography (LLC), the dominant interaction field is created by the solvent system composition itself. The system is therefore tuned by choosing the two-phase composition, hydrophobicity, pH, ionic strength, additives, ligands or dissolved selectors and by defining which phase is stationary and which is the mobile phase. Because phase roles can be inverted, a single solvent system can often generate mirrored normal-phase-like and reversed-phase-like behavior, depending on the application or the economy-driven target of the final product. Changing between operation mode (ascending or descending) can influence the final solvent next to the Col, the elution order of critical impurities vs. Col to minimize cross contamination.

In addition to the operational mode options, LLC can also offer further setup similar to MCSGP (multicolumn countercurrent solvent gradient purification), SMB (simulated moving bed) or other continuous chromatographic operation as well [7–9].

The liquid–liquid platform does not reproduce every liquid–solid mode one-to-one, but it does address a surprisingly large fraction of the same industrial problem space through different physical levers, each with a distinct industrial application [18–21]. More broadly, liquid–liquid partition chromatography has also been extended into continuous and hybridized process configurations, reinforcing that its relevance is not limited to classical batch operation [18]. At the same time, any implementation of CPC in place of preparative HPLC should be understood not as a simple equipment substitution, but as a manufacturing process change requiring comparability, validation and an appropriate control strategy [19–21]. pH-zone refining is a displacement-like technique in which a titratable sample is loaded at high concentration and separated by pH gradient across the two phases; it is particularly powerful for ionizable compounds at preparative and semi-industrial scale, where it can produce sharp, concentrated fractions at high loading [22, 23]. Elution-extrusion addresses the practical problem of strongly retained components: once elution is complete, the stationary phase itself is pushed out, recovering material that would otherwise remain bound in a conventional column [7, 9]. Increasing the applicable polarity range when compared to traditional column chromatography gradient setup. Furthermore, extraction protocol at the end of each loop limits potential batch-to-batch carryover and provides consistent performance throughout the entire process since the immobilized liquid acting as stationary phase is replenished between each run, eliminating the slow decay of stationary phase performance, experienced in HPLC regime.

Dual-mode protocols exploit phase inversion to change selectivity mid-run, allowing the same solvent system to separate compounds with very different partition coefficients in a single operation [7]. Charge-sensitive separations not amenable to standard pH control can be addressed through ion-pairing or ionic selectors dissolved directly in one phase, without need for ion-exchange media. Chiral discrimination and affinity-like selectivity can both be introduced through selectors or ligands dissolved in one phase of the biphasic system, rather than through immobilized functionalities bound to a solid matrix, which means the chiral selector loading and identity can be varied between runs without column replacement. Aqueous two-phase systems (ATPS) extend the platform into gentle biomolecule-oriented partitioning space, where polyethylene glycol–salt or polymer–polymer systems can selectively partition proteins, nucleic acids, or cell organelles without organic solvent exposure [9].

It is equally important to note the practical limitation common to all CPC solvent system development: unlike HPLC, where extensive stationary-phase libraries, retention prediction tools, and method databases exist, solvent system selection in CPC remains substantially more empirical.

To address this technical gap, RotaChrom and NLR (National Laboratories of the Rockies, US) are developing solvent libraries and computational workflows to support CPC solvent system selection. NRL has reported an HPC-enabled (high performance computing) tool for *in silico* high-throughput screening of CPC solvent systems, while recent NLR work has combined partition coefficient measurements with CUP (cell utilized partitioning) model-based refinement to guide solvent-system selection and optimization [24–26].

3. From retention mechanism to purification duty

Industrial purification in sequence must tolerate dirty feeds, high impurity burden, low solubility margins, and high mass throughput. Afterwards, discriminating between structurally similar species at low relative abundance while preserving final purity, recovery and control. Two further process metrics deserve explicit attention here, because they govern industrial platform selection as much as separation performance does: loading capacity and cycle productivity. In preparative HPLC, column loading is constrained by competitive displacement and stationary phase saturation; at preparative and pilot scale, typical mass loadings in overloaded gradient mode depending on the system and solubility in the mobile phase often imposes an additional practical ceiling.

CPC loading is governed by the partition coefficient and the volume of retained stationary phase and because there is no solid surface to saturate, the system can in principle tolerate crude, concentrated feeds without the peak shape distortion associated with solid-phase overloading. At preparative scale, reported CPC mass loadings for favorable systems commonly reach higher concentration in stationary phase [7, 22] and loading tolerance for complex matrices is generally broader than in packed-bed systems. The scalability principle of CPC instruments is engineered for linear scale-up, meaning that separation performance, retention and selectivity developed at laboratory scale translate predictably to larger rotors by maintaining the same S_r , flow rate per unit volume and rotational speed [7, 9]. This is a significant industrial advantage. In preparative HPLC, scale-up is not linear: bed geometry, packing uniformity, pressure drop and wall effects all change non-trivially as column diameter increases and substantial re-optimization is typically required at each scale transition.

In CPC, the rotor volume is the primary scale variable and the separation logic does not change with it. This means that a method developed at 250 mL rotor scale can be transferred to industrial scale instruments with high confidence in outcome, substantially reducing scale-up risk and development time. Productivity, expressed as mass of purified product per unit time and system volume, is a more complex comparison. CPC flow rates are constrained by the need to retain the stationary phase under centrifugal force and are generally lower per unit volume than those achievable in modern preparative HPLC.

For resolution-demanding, low-selectivity separations, this flow rate limitation reinforces the efficiency gap and makes HPLC the higher-throughput choice, though the productivity can only be enhanced if the bed size increases significantly as the concentration of the feed kept at low level (which is directly impacting CAPEX/OPEX). For high-loading, high-selectivity separations where CPC's partition-driven discrimination is sufficient, the productivity comparison favors CPC, particularly when HPLC cycle time includes solvent switching, re-equilibration and column regeneration steps. Neither platform dominates productivity universally, the comparison depends on the specific selectivity, loading and operational context of each separation (Fig. 4.).




Figure 4. Decision framework for choosing CPC over preparative HPLC. Selection branches on solvent-system selectivity, then on loading/surface-sensitivity or resolution demand. High-loading, high-selectivity or adsorption-sensitive separations favor CPC; resolution-demanding, low-selectivity separations favor HPLC, though its throughput gain is conditional on larger bed volumes at low feed concentration (CAPEX/OPEX impact).

The liquid–solid route is often attractive when high efficiency, mature method culture and a clearly defined stationary phase chemistry are required. The liquid–liquid route is often attractive when the purification problem benefits from solvent system flexibility, thus selectivity engineering, orthogonal selectivity, broad loading tolerance, or avoidance of uncontrolled interactions with solid surfaces. Neither route is inherently confined to a supporting role, both can function as stand-alone purification systems when their governing mechanism and separation logic matches the impurity problem and process context.

For that reason, comparison should not begin and end with nominal plate number. A higher theoretical plate number is valuable, but it does not by itself eliminate surface-driven losses, loading constraints, solubility mismatch, pressure-drop burden, or scale-up fragility. The real question is whether the separation problem is limited by insufficient selectivity or by the cost of achieving resolution. If α remains too low, there is no real substitute for HPLC: the separation depends on high efficiency and CPC is unlikely to compensate for the gap in N . If α is sufficiently large, however, the logic changes. In that regime, CPC may outperform HPLC not by generating more plate number but by delivering the required purification with lower burden, broader loading capacity and more favorable CAPEX/OPEX economics.

4. Stationary phase burden in packed bed preparative systems

In packed-bed chromatography, the solid matrix is not a passive support. Its particle morphology, pore architecture, surface chemistry, packing behavior, chemical stability and lifetime under use all become process variables. At analytical scale, these constraints may remain manageable. At process scale, they translate directly into pressure, downtime, cleaning restrictions, replacement burden, reproducibility risk and cost. The comparison between HPLC and CPC is whether the separation depends on the continuous packing, protection, cleaning, regeneration and eventual replacement of a solid stationary phase at all [1, 3–6, 12, 13].

In the following sections the main pitfalls of the stationary phase are detailed, which reinforces the need for careful evaluation of a proper purification strategy.


4.1. Silica: high efficiency built on narrow chemical and manufacturing tolerances

Silica remains the dominant matrix of HPLC because it offers exactly the combination that makes packed-bed efficiency possible: rigid particles, high surface area, and chemically modifiable surfaces. Yet this performance rests on a narrow and demanding manufacturing window. High-efficiency operation requires highly spherical particles with tight size distribution, but producing such material reproducibly is expensive, especially as particle size decreases. Smaller particles improve plate count, but they also intensify manufacturing difficulty, increase pressure requirements during use and raise the likelihood of fines generation. Those fines are not merely a production issue; later they contribute to frit clogging, disturbed flow and pressure rise inside the column [4, 12, 13].

A second limitation lies in pore architecture. Silica performance depends on balancing pore size, pore uniformity and internal surface area, but these variables do not scale independently. High surface area is typically associated with smaller pores, whereas larger pores required for larger molecules reduce accessible area. Maintaining a consistent pore diameter across a batch is technically demanding and internal functionalization becomes harder as pores narrow and access to the inner surface becomes more restricted. Silica therefore always operates inside a structural compromise: it can be optimized but not maximized in every direction at once [4, 5].

Chemical stability imposes a further constraint. Standard silica is not broadly stable across the full pH range relevant to many industrial purification environments. Under alkaline conditions the matrix gradually dissolves, whereas under strongly acidic conditions bonded phases may degrade or cleave. Residual silanol groups remain another persistent limitation. Even after end-capping, they contribute to secondary interactions, peak tailing, nonspecific adsorption and progressive fouling. Moisture control adds yet another layer of sensitivity, since silica surface hydration influences chromatographic behavior and must be controlled tightly if reproducibility is maintained.

At analytical scale, those tolerances are often acceptable; at process scale, they become operational constraints [4, 5].



4.2. Process-scale HPLC: when the packed bed becomes the burden

Despite its high resolving power and established industrial status, preparative HPLC becomes increasingly burdened as throughput and scale are increased. Large-diameter packed beds must be produced and maintained with high uniformity if acceptable efficiency and asymmetry are to be preserved. For this reason, compression-based column designs such as dynamic axial compression (DAC) were introduced to stabilize the bed and compensate for packing relaxation during operation. Even so, bed packing remains a critical variable: packing non-uniformity promotes peak asymmetry, void formation may develop during continued operation and the development of robust packing procedures at large scale is time-consuming and costly. The bed therefore ceases to be just a medium for selectivity and becomes a mechanical system that must be packed reproducibly, stabilized under flow, and protected against degradation. This is the first major burden of packed-bed preparative chromatography: the separation may still be chemically powerful, but the bed itself must continuously be defended against mechanical failure [12, 13].

Additional limitations arise from feed handling and repeated column use. In preparative RP-HPLC, limited solubility in mobile-phase-compatible solvents often forces the use of large injection volumes, while volume overload may distort peak shape and reduce effective utilization of the bed [14]. This is especially problematic when crude feeds are impurity-rich and purification must be repeated across many cycles. What begins as a loading problem quickly becomes a throughput and process-economics problem [3, 14].

At the same time, column performance cannot be assumed to remain constant indefinitely. Regeneration, cleaning, re-equilibration and in some cases repacking or media replacement are required if acceptable chromatographic performance is to be maintained [3, 12, 13]. HPLC therefore remains a powerful high-resolution platform, but at industrial scale its packed-bed architecture imposes a significant burden through maintenance, downtime and process economics.

Moreover, at industrial scale, the burden of preparative HPLC extends beyond separation performance to substantial solvent consumption and waste generation. RP-HPLC is commonly operated with aqueous buffers and organic modifiers, most often acetonitrile and solvent demand increases markedly with scale, repeated cycling, washing, re-equilibration and column cleaning. Solvent-rich waste streams therefore add both economic and environmental burden through storage, recovery, segregation, treatment and disposal; and have been identified as major contributors to process mass intensity and downstream sustainability challenges [15–17].

Taken together, these limitations do not diminish the value of HPLC as a high-resolution industrial platform, but they do clarify where its logic becomes burdensome: the separation performance of the method is inseparable from a packed solid bed that must be packed, stabilized, protected, cleaned, regenerated and eventually replaced. At scale, the question is therefore no longer only whether HPLC can deliver the required purity, but at what process cost that purity is achieved [3, 12–17].

4.3. Polymeric beads and monoliths: different matrices, different liabilities

Silica is not the only stationary phase option at process scale, but the alternatives do not remove the trade-off; they shift it. Polymeric beads dominate bioprocessing because they solve problems that silica handles poorly. Their broader pH tolerance allows aggressive clean-in-place procedures that silica cannot survive. Their hydrated three-dimensional structure gives large biomolecules more usable internal binding volume and their hydrophilic background reduces nonspecific adsorption and lowers the risk of denaturing fragile proteins. In biologics manufacturing, these are decisive advantages.

However, polymeric matrices solve one burden while intensifying another. Their principal weakness is mechanical softness. Under higher flow and pressure, many polymeric beds deform, compact or shift. Interstitial space collapses, pressure rises and the bed becomes less stable under precisely the conditions where higher throughput would be desirable. Their separation efficiency is also typically lower than that of silica because polymeric beads deliver broader peaks and fewer theoretical plates. Mass transfer remains diffusion-limited, which often forces slower flow rates if acceptable separation is to be maintained. Synthetic polymers introduce an additional instability through solvent-dependent swelling and shrinking, while natural polymer matrices add storage and microbial vulnerabilities that must also be managed [1–3].

Monoliths reformulate the problem again. Because they replace packed particles with a continuous porous structure, they offer convective mass transport and low pressure drop. This makes them attractive for very large biomolecules such as plasmids, viral vectors and related systems where diffusion-limited bead beds struggle. They also avoid classical slurry packing, which is a real operational advantage. Yet monoliths carry their own constraints. Large-scale manufacture is difficult because *in situ* polymerization generates exothermic heat that is hard to dissipate evenly, leading to structural defects, pore heterogeneity and cracking. For peptides and many smaller or mid-sized targets, monoliths often lack sufficient binding capacity because their usable surface area is too limited. They are also unforgiving to fouling: once the front face is blocked, there is usually no practical unpacking or partial remediation route. Instead of packing burden, the process inherits a disposable-module cost burden [1, 2].

4.4. Performance decay: why packed-bed performance does not remain static

A final limitation is reproducibility, which is in reality a lifecycle issue: process-scale packed-bed chromatography does not operate at fixed performance. From the first crude injection onward, the column enters a state of progressive drift.

Retention times shift because strongly retained impurities progressively blind active sites. Selectivity narrows because fouled regions no longer behave like the fresh bed originally validated. Peak shapes broaden and tail as frit fouling, pore blockage, ligand loss and bed disturbance degrade mass-transfer behavior. Baseline pressure rises because fines accumulate or the bed physically deforms. Yield then becomes entangled with purity control: to maintain specification, collection windows must be tightened and as they are tightened, recovery falls. If they are not tightened, purity drifts instead [3, 12, 13].

This means that a process-scale HPLC column often becomes economically worse before it becomes technically unusable. It may still function chemically, but no longer efficiently enough to justify the solvent use, cycle count, labor and recovery loss required to keep the process within specification. Performance decay is therefore not just an operational inconvenience. It is part of process economics [3, 12, 13, 15–17].

Manufacturing manages this through lifecycle qualification rather than any assumption of static performance. Columns are validated conservatively, monitored through pressure, plate-height and asymmetry trends and retired before catastrophic failure. But that does not remove the burden; it formalizes it. Every repacking event means unpacking, cleaning, slurry preparation, repacking, validation and the risk of failed qualification. At large scale, this is measured in days of disruption rather than hours [12, 13].

At scale, the question is therefore no longer only whether HPLC can deliver the required purity, but at what process cost that purity is achieved.

5. Conclusion

CPC and HPLC are not interchangeable tools, but different separation logics. The real industrial question is that which one is carrying the right purification duty. In many downstream trains, that means keeping HPLC where very high resolution is genuinely needed, while using CPC where matrix-free, partition-driven separation can remove burden earlier, tolerate higher loading and improve the overall economics of purification.

Contact

Interested in the method, achievable purities, or fraction portfolio? Contact us to discuss your target alkaloids, required specifications and scale.

ROTACHROM.COM
ROTACHROM.COM/CONTACT

RotaChrom's CPC technology delivers highly selective separations making it well suited for complex compounds with similar physical and chemical properties across pharma, biotech, botanicals, and nutraceuticals. It fits naturally into downstream processing workflows, bridging the gap between bulk extraction and high-purity isolation where extraction alone is no longer sufficient to meet purity or regulatory requirements.

From APIs and bioproducts to natural bioactives and functional ingredients, CPC ensures high yield, reproducibility and gentle handling of sensitive molecules.

RotaChrom's technology is built for true scalability, preserving the same separation principles from benchtop development through pilot-scale to full industrial production. The same methodology, solvent systems and process logic can be transferred without re-engineering, minimizing scale-up risk and validation effort.



References

- [1] Guiochon G. Preparative liquid chromatography. *J Chromatogr A*. **2002**; 965(1–2):129–161. doi:10.1016/S0021-9673(01)01471-6.
- [2] Rajendran A. Recent developments in preparative chromatographic processes. *Curr Opin Chem Eng*. **2013**; 2(2):263–270. doi:10.1016/j.coche.2013.04.001.
- [3] Andersson L, Persson J. *Aspects of industrial purification of peptides using large-scale chromatography*. Limhamn (Sweden): PolyPeptide Laboratories (Sweden) AB; 2000.
- [4] Qiu H, Liang X, Sun M, Jiang S. Development of silica-based stationary phases for high-performance liquid chromatography. *Anal Bioanal Chem*. **2011**; 399(10):3307–3322. doi:10.1007/s00216-010-4611-x.
- [5] Bocian S, Vajda P, Felinger A, Buszewski B. Residual silanols at reversed-phase silica in HPLC: a contribution for a better understanding. *J Sep Sci*. **2012**; 35(10–11):1191–1200. doi:10.1002/jssc.201200055.
- [6] Guimarães GJ, Bartlett MG, et al. Managing nonspecific adsorption to liquid chromatography hardware: a review. *J Sep Sci*. **2023**. doi:10.1002/jssc.202300060.
- [7] Bojczuk M, Nowak D, Olech M. Centrifugal partition chromatography: a review of recent applications and some classic references. *J Sep Sci*. **2017**; 40(7):1597–1609. doi:10.1002/jssc.201601221.
- [8] Liu Y, Friesen JB, McAlpine JB, Pauli GF. Solvent system selection strategies in countercurrent separation. *Planta Med*. **2015**; 81(18):1582–1591. doi:10.1055/s-0035-1558110.
- [9] Berthod A. Separation and purification with a liquid stationary phase. *Separations*. **2017**; 4(4):30. doi:10.3390/separations4040030.
- [10] Ettre LS, Csaba Horváth The development of the first modern high-performance liquid chromatograph. *LCGC North Am*. **2005**; 23(5).
- [11] Yale News. In memoriam: pioneering leader who developed widely used tools of chemical analysis, Csaba Horváth. New Haven (CT): Yale University; 2004 Apr 19.
- [12] Williams JA, Przybycien TM, Cheng YL. Maintenance of column performance at scale. *J Chromatogr A*. **2002**; 971(1–2):211–221.
- [13] Dorn M, Eschbach F, Hekmat D, Weuster-Botz D. Influence of different packing methods on the hydrodynamic stability of chromatography columns. *J Chromatogr A*. **2017**; 1516:89–101. doi:10.1016/j.chroma.2017.08.019.
- [14] Miller L, Ayube SJ. In-column dilution: improving volume loadability in reverse-phase gradient chromatography through the use of a silica pre-column. *J Chromatogr A*. **2020**; 1618:460897. doi:10.1016/j.chroma.2020.460897.
- [15] Yabré M, Ferey L, Somé IT, Gaudin K. Greening reversed-phase liquid chromatography methods using alternative solvents for pharmaceutical analysis. *Molecules*. **2018**; 23(5):1065. doi:10.3390/molecules23051065.
- [16] Passerini BFG, Breitzkreitz MC. Recent developments in green chromatography. *Braz J Anal Chem*. **2023**; 10:17–34. doi:10.30744/brjac.2179-3425.RV-126-2022.
- [17] Shen Y, Chen B, van Beek TA. Alternative solvents can make preparative liquid chromatography greener. *Green Chem*. **2015**; 17:4073–4081. doi:10.1039/C5GC00887E.
- [18] Viano R, Milani L, Örkényi R, Adamo A. Continuous separation and purification of chemical species by continuous liquid-liquid partition chromatography implemented with membrane-based phase separators. *Org Process Res Dev*. **2023**; 27(12):2296–2308. doi:10.1021/acs.oprd.3c00239.
- [19] International Council for Harmonisation (ICH). *ICH Q8(R2): Pharmaceutical development*. Geneva: ICH; 2009.
- [20] U.S. Food and Drug Administration. *PAT – A framework for innovative pharmaceutical development, manufacturing, and quality assurance: guidance for industry*. Silver Spring (MD): FDA; 2004.
- [21] U.S. Food and Drug Administration. *Process validation: general principles and practices: guidance for industry*. Silver Spring (MD): FDA; 2011.
- [22] Berthod A, Maryutina T, Spivakov B, Shpigun O, Sutherland IA. Countercurrent chromatography in analytical chemistry (IUPAC Technical Report). *Pure Appl Chem*. **2009**; 81(2):355–387. doi:10.1351/PAC-REP-08-06-05.
- [23] Ito Y. pH-zone-refining countercurrent chromatography: a preparative and analytical tool. *J Chromatogr A*. **1996**; 753(1):1–36. doi:10.1016/S0021-9673(96)00300-5.
- [24] National Renewable Energy Laboratory (NREL). *Advanced Computing Annual Report 2021*. Golden (CO): NREL; 2022. Report No.: NREL/PR-2C00-82265.
- [25] Choi H, Alherech M, Jang JH, Woodworth SP, Ramirez KJ, Karp EM, Beckham GT. Counter-current chromatography for lignin monomer–monomer and monomer–oligomer separations from reductive catalytic fractionation oil. *Green Chem*. **2024**; 26:5900–5914. doi:10.1039/D4GC00765D.
- [26] RotaChrom Technologies PLC. *Experimental design* [Internet]. Kecskemét: RotaChrom Technologies PLC; [cited 2026 May 29]. Available from: https://rotachrom.com/services-solution_planning_experimental_design/